**Technical Specification**

Part No.	Alloy	Adapter Style	Product Weight	Build Height	Max Patient Weight	Screw Torque
4T08	Titanium	Receiver	48 g	12 mm.	136 Kg. / 300 lbs.	15 Nm.
4T09	Titanium	Pyramid	38 g	8 mm.	136 Kg. / 300 lbs.	
4T10	Titanium	Pyramid	44 g	10 mm.	136 Kg. / 300 lbs.	
4T11	Titanium	Pyramid	48 g	12 mm.	136 Kg. / 300 lbs.	
3T90	Titanium	Frame	33 g	9 mm.	136 Kg. / 300 lbs.	9 Nm.
4S35	S. Steel	Frame	109 g	9 mm.	136 Kg. / 300 lbs.	9 Nm.
4A80	Aluminum	Plate	70 g	10 mm.	136 Kg. / 300 lbs.	

**Description and Purpose**

These instructions are for use by the practitioner. This product is designed to be used for lower limb external prosthetic length, toe direction and alignment corrections. Intended for single user only. This product can be used for mobility classes 1-4. There are no known contra-indications if used in accordance with these instructions. Ensure that the user has understood any instructions for use, drawing particular attention to the maintenance information.

**Construction**

Use 4mm Allen wrench for all screws. Before fitting this product onto the patient, all parts must be joined properly and any set screws must be secured enough for the patient to walk. In order to make final adjustments, there should be no loose screws. After the final adjustments, apply single drop of removable Loctite® 242 to each of the set/cap screws and then torque the product.

**Function**

Use 4mm Allen wrench for all clamp screws. Tape the open areas of the 3 arm anchor frame (4S35) or use silicone putty to keep lamination away from threads. 4 hole frame (3T90) connect directly to 4A80 lamination plate with M6 countersunk attachment screws (12Nm). 6mm mounting screw closest to the 5mm clamp screw (on 3T90) must be left loose until the rotation is set. Once the rotation set screw is tightened to torque specification listed above, go back and tighten the last mounting screw. Improper installation or uneven, insufficient or excessive torque may cause bolt or adapter failure. Always use quality metric tools for installation. When tightened securely using Loctite on the pyramid screws only and torque settings shown, the device offers a secure way to connect two parts of the limb.

**Limitations on Use**

Intended life:  
 A local risk assessment should be carried out based upon activity and usage.

Lifting loads:  
 User weight and activity is governed by the stated limits. Load carrying by the user should be based on a local risk assessment.

Environment:  
 Thoroughly rinse with fresh water after use in abrasive environments such as those that may contain sand or grit, for example, to prevent wear or damage to moving parts. Thoroughly rinse with fresh water after use in salt or chlorinated water. Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).

**Maintenance**

Maintenance is to be carried out by competent personnel only. An annual visual inspection is recommended. During maintenance, it is appropriate to consider any weight or activity increase by the user. If the prosthesis is exposed to circumstances or impact for which it was not designed, inspection by a qualified practitioner is recommended.

- Ensure all screws are secure. If not, remove and clean screws, reapply Loctite to each of the set/cap screws, tighten all screws to the correct torque settings.
- Check for visual defects that may affect proper function.

The user should be advised to contact their practitioner if:

- Their condition changes.
- Performance of the device changes.
- Changes in performance may include:
  - Any unusual noises
  - Excessive play or loss of alignment

**Warranty**

All MEDEX Modular Adapters are unconditionally guaranteed for 2 years. Cutting, drilling, modifying or using the component other than as described in this manual will void the warranty.

**Liability**

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device.

**CE Conformity**

This product meets the requirements of MDR 2017/745 guidelines for medical products. This product has been classified as a Class 1 Product according to the classification criteria outlined in Appendix IX of the guidelines.

MEDEX INTERNATIONAL, INC. – 2720 Duvall Road, Burtonsville, MD 20866 – Phone: 1-888-886-2420 – Fax: 301-657-2796 – Last Revised 26/5/21

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Artikelnr.	Material	Adaptertyp	Gewicht	Bauhöhe	Gewichtslimit	Anzugsmoment
4T08	Titan	Aufnahme	48 g	12 mm.	136 Kg. / 300 lbs.	15 Nm.
4T09	Titan	Pyramide	38 g	8 mm.	136 Kg. / 300 lbs.	
4T10	Titan	Pyramide	44 g	10 mm.	136 Kg. / 300 lbs.	
4T11	Titan	Pyramide	48 g	12 mm.	136 Kg. / 300 lbs.	
3T90	Titan	Rahmen	33 g	9 mm.	136 Kg. / 300 lbs.	9 Nm.
4S35	Edelstahl	Rahmen	109 g	9 mm.	136 Kg. / 300 lbs.	9 Nm.
4A80	Aluminium	Scheibe	70 g	10 mm.	136 Kg. / 300 lbs.	

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